

Inspired by Nature.
Perfect by Design.

ZirLife™ AESTHETIC


TRANSLUCENT ZIRCONIA CROWN

CHAIRSIDE TABLE GUIDE


- Indicated for anterior crowns & bridges
- 803 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Dipped stained for true shades
- Polished after milling to remove micro ridges



©2019 CLASSIC CRAFT DENTAL LABORATORY



ANTERIOR CROWNS



ANTERIOR BRIDGES

SEE POSTERIOR PREP ON OTHER SIDE

Classic★Craft
DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007
2807 Roger Lacey Drive #B, Savannah, GA 31404
classiccraftdental.com

