

Inspired by Nature.
Perfect by Design.

ZirLifeUSA™

FULL CONTOUR ZIRCONIA CROWN

CHAIRSIDE TABLE GUIDE

- Indicated for posterior crowns & bridges
- >1100 mpa strength
- Medium Translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Dipped stained for true shades
- Polished after milling to remove micro ridges




©2019 CLASSIC CRAFT DENTAL LABORATORY

Classic★Craft

DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007
2807 Roger Lacey Drive #B, Savannah, GA 31404
classiccraftdental.com



POSTERIOR CROWNS

SEE ANTERIOR PREP ON OTHER SIDE

