

# Scan bodies, model analogs and printed models – Compatibility chart



Manufacturer	Implant	Atlantis FLO	Atlantis IO FLO	Elos Accurate IO Scan Body Only Atlantis Implant Kit for 3Shape (DME)	Elos Accurate Model Analog Only Atlantis Implant Kit for 3Shape (DME)	Atlantis Printed Model
		Code (REF)	Code (REF)	REF	REF	Supported
<b>BioHorizons</b>	Internal, Tapered Internal 3.5	B-01 (34527)	IO-B-01 (35212)	IO 5A-A	MA-ZSV35-1	YES
	Internal, Tapered Internal 4.5	B-02 (34528)	IO-B-02 (35213)	IO 5A-B	MA-ZSV45-1	YES
	Internal, Tapered Internal 5.7	B-03 (34529)	IO-B-03 (35214)	IO 5A-C	MA-ZSV57-1	YES
	External Implant 3.5	A-03 (34518)	–	–	–	–
	External Implant 4.0	A-06 (34521)	–	–	–	–
	External Implant 5.0	A-09 (34524)	–	–	–	–
<b>Biomet 3i</b>	Certain 3.25, 4/3	F-01 (34539)	IO-F-01 (35225)	IO 7B-A	MA-BCE34-1	YES
	Certain 4.0, 5/4	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	MA-BCE41-1	YES
	Certain 5.0	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	MA-BCE50-1	YES
	Certain 6.0	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	MA-BCE60-1	YES
	Certain Prevail 3/4/3 (3.4), 4/3	F-01 (34539)	IO-F-01 (35225)	IO 7B-A	MA-BCE34-1	YES
	Certain Prevail 4/5/4, 5/4	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	MA-BCE41-1	YES
	Certain Prevail 5/6/5, 6/5	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	MA-BCE50-1	YES
	Certain XP 4/5	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	MA-BCE50-1	YES
	Certain XP 5/6	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	MA-BCE60-1	YES
	Miniplant 3.25	A-07 (34522)	–	–	–	–
	Osseotite 3.75, 4.0	A-07 (34522)	–	–	–	–
	Osseotite 5.0, 6.0	A-07 (34522)	–	–	–	–
	Osseotite XP 4/5	A-07 (34522)	–	–	–	–
	Osseotite XP 5/6	A-07 (34522)	–	–	–	–

Please turn over →

Atlantis FLO, Atlantis IO FLO and Elos Accurate IO Scan Body can be used with all connectivity when scanning for Atlantis abutments and are also included in the Atlantis Implant Kit for 3Shape (DME). 3Shape Scan Bodies and Elos Accurate Model Analogs (Gen 1) are also included in the DME file. This chart is valid for DME file created on 05-14-2019 or later.

**Note:** 3Shape Scan Bodies can only be used when sending in the scan files using the validated direct workflow from 3Shape to Atlantis WebOrder. The scan body detection must be done in the 3Shape software, therefore it is **NOT** possible to scan with 3Shape Scan Bodies, export the .stl files and send the files using Atlantis Scan Upload or any other .stl upload.

Atlantis Printed Model is available to order in Atlantis WebOrder and includes the Elos Accurate Model Analog Gen 2.

All trademarks and company names are the property of their respective owners.

Manufacturer	Implant	Atlantis FLO	Atlantis IO FLO	Elos Accurate IO Scan Body Only Atlantis Implant Kit for 3Shape (DME)	Elos Accurate Model Analog Only Atlantis Implant Kit for 3Shape (DME)	Atlantis Printed Model
		Code (REF)	Code (REF)	REF	REF	Supported
<b>Camlog</b>	CONELOG 3.3	Q-01 (34666)	—	—	—	—
	CONELOG 3.8	Q-02 (34667)	—	—	—	—
	CONELOG 4.3	Q-02 (34667)	—	—	—	—
	CONELOG 5.0	Q-03 (34668)	—	—	—	—
	Screw-line Implant 3.3	N-01 (34584)	IO-N-01 (35239)	IO 9A-A	MA-CAS33-1	YES
	Screw-line/Root-line Implant 3.8	N-02 (34585)	IO-N-02 (35240)	IO 9A-B	MA-CAS38-1	YES
	Screw-line/Root-line Implant 4.3	N-03 (34586)	IO-N-03 (35241)	IO 9A-C	MA-CAS43-1	YES
	Screw-line/Root-line Implant 5.0	N-04 (34587)	IO-N-04 (35242)	IO 9A-D	MA-CAS50-1	YES
	Screw-line/Root-line Implant 6.0	N-04 (34587)	IO-N-04 (35242)	IO 9A-D	—	YES
<b>Dentsply Sirona</b> <small>Note: OsseoSpeed implants are part of the Astra Tech Implant System.</small>	Ankylos C/X 3.5, 4.5, 5.5, 7.0 (non-indexed)	L-01 (34583)	IO-L-01 (35235)*	IO 8B-A SA	MA-DAN35/70-1	YES
	OsseoSpeed EV 3.0  GREEN	P-01 (34593)	IO-P-01 (35243)	IO 3B-A SA	MA-ATE30-1**	YES
	OsseoSpeed EV 3.6  PURPLE	P-02 (34594)	IO-P-02 (35244)	IO 3B-B SA	MA-ATE36-1**	YES
	OsseoSpeed EV 4.2  YELLOW	P-03 (34595)	IO-P-03 (35245)	IO 3B-C SA	MA-ATE42-1**	YES
	OsseoSpeed EV 4.8  BLUE	P-04 (34596)	IO-P-04 (35246)	IO 3B-D SA	MA-ATE48-1**	YES
	OsseoSpeed EV 5.4  BROWN	P-05 (34597)	IO-P-05 (35247)	IO 3B-E SA	MA-ATE54-1**	YES
	OsseoSpeed Profile EV 4.2  YELLOW	P-06 (34598)	IO-P-06 (35248)	IO 3C-A SA	—***	YES
	OsseoSpeed Profile EV 4.8  BLUE	P-07 (34599)	IO-P-07 (35249)	IO 3C-B SA	—***	YES
	OsseoSpeed / OsseoSpeed TX 3.0  YELLOW	C-01 (34530)	IO-C-01 (35215)	IO 3A-A SA	MA-ATO30-1	YES
	OsseoSpeed / OsseoSpeed TX 3.5, 4.0  AQUA	C-02 (34531)	IO-C-02 (35216)	IO 3A-B SA	MA-ATO35-1	YES
	OsseoSpeed / OsseoSpeed TX 4.5, 5.0  LILAC	C-03 (34532)	IO-C-03 (35217)	IO 3A-C SA	MA-ATO45-1	YES
	OsseoSpeed TX Profile 4.5, 5.0	C-04 (34592)	IO-C-04 (35218)	—	—	—
	Xive S 3.0  3.0	J-01 (34575)	IO-J-01 (35230)	IO 8A-A	MA-DX130-1	YES
	Xive S 3.4  3.4	J-02 (34576)	IO-J-02 (35231)	IO 8A-B	MA-DX134-1	YES
	Xive S 3.8  3.8	J-03 (34577)	IO-J-03 (35232)	IO 8A-C	MA-DX138-1	YES
	Xive S 4.5  4.5	J-03 (34577)	IO-J-04 (35264)	IO 8A-D	MA-DX145-1	YES
Xive S 5.5  5.5	J-03 (34577)	IO-J-05 (35265)	IO 8A-D	—	YES	
<b>Hiossen</b>	ET Mini 3.5	R-01 (34671)	—	—	—	—
	ET Regular 4.0, 4.5, 5.0, 6.0, 7.0	R-02 (34672)	—	—	—	—
<b>Keystone Dental</b>	PrimaConnex SD 3.3, Genesis 3.8	M-01 (34589)	—	—	—	—
	PrimaConnex SD 3.5	M-01 (34589)	—	—	—	—
	PrimaConnex RD 4.0, Genesis 4.5	M-02 (34590)	—	—	—	—
	PrimaConnex RD 4.1	M-02 (34590)	—	—	—	—
	PrimaConnex WD 5.0, Genesis 5.5, 6.5	M-03 (34591)	—	—	—	—

Please turn over

All trademarks and company names are the property of their respective owners.

\* The Atlantis IO FLO for Ankylos C/X is for single use only.

\*\* Astra Tech Implant System Lab Abutment Screw EV is designed for use with Astra Tech Implant System Implant Replica EV only and will not fit into the ELOS Accurate Model analog. Either order an extra Atlantis abutment screw, Abutment Screw EV or use ELOS lab screw.

\*\*\* Scan body/Model analog exist but is not supported for Atlantis solutions.

Manufacturer	Implant	Atlantis FLO	Atlantis IO FLO	Elos Accurate IO Scan Body Only Atlantis Implant Kit for 3Shape (DME)	Elos Accurate Model Analog Only Atlantis Implant Kit for 3Shape (DME)	Atlantis Printed Model
		Code (REF)	Code (REF)	REF	REF	Supported
<b>MIS</b>	M4 NP 3.3	B-04 (34678)	—	—	—	—
	M4 SP 3.75, 4.2	B-01 (34527)	IO-B-01 (35212)	—	—	—
	M4 WP 5.0, 6.0	B-02 (34528)	IO-B-02 (35213)	—	—	—
	SEVEN NP 3.3	B-04 (34678)	—	—	—	—
	SEVEN SP 3.75, 4.2	B-01 (34527)	IO-B-01 (35212)	—	—	—
	SEVEN WP 5.0, 6.0	B-02 (34528)	IO-B-02 (35213)	—	—	—
<b>Nobel Biocare</b>	NobelActive 3.0	K-03 (34665)	—	IO 2B-D SA	MA-NBA30-1	YES
	NobelActive / NobelReplace Conical Connection NP 3.5	K-01 (34578)	IO-K-01 (35233)	IO 2B-A SA	MA-NBA35-1	YES
	NobelActive / NobelReplace Conical Connection RP 4.3	K-02 (34579)	IO-K-02 (35234)	IO 2B-B SA	MA-NBA43-1	YES
	NobelActive / NobelReplace Conical Connection RP 5.0	K-02 (34579)	IO-K-02 (35234)	IO 2B-B SA	MA-NBA50-1	YES
	NobelActive WP 5.5	K-04 (34677)	—	—***	—***	—
	NobelReplace NP 3.5	D-01 (34533)	IO-D-01 (35219)	IO 2A-A	MA-NBR35-1	YES
	NobelReplace RP 4.3	D-02 (34534)	IO-D-02 (35220)	IO 2A-B	MA-NBR43-1	YES
	NobelReplace WP 5.0	D-03 (34535)	IO-D-03 (35221)	IO 2A-C	MA-NBR50-1	YES
	NobelReplace 6.0	D-04 (34536)	IO-D-04 (35222)	IO 2A-D	MA-NBR60-1	YES
	Brånemark System NP 3.3	A-04 (34519)	—	IO 6A-A	MA-BRA35-1	YES
	Brånemark System RP 3.75, 4.0	A-06 (34521)	—	IO 6A-B	MA-BRA41-1	YES
Brånemark System WP 5.0, 5.5	A-11 (34526)	—	IO 6A-C	MA-BRA51-1	YES	
<b>Straumann</b>	Bone Level 3.3 NC (Narrow CrossFit)	H-01 (34542)	IO-H-01 (35228)	IO 4A-A SA	MA-SBO33-1	YES
	Bone Level 4.1 RC (Regular CrossFit)	H-02 (34543)	IO-H-02 (35229)	IO 4A-B SA	MA-SBO41-1	YES
	Bone Level 4.8 RC (Regular CrossFit)	H-02 (34543)	IO-H-02 (35229)	IO 4A-B SA	MA-SBO48-1	YES
	Standard/Standard Plus 4.8 RN (Regular Neck)	E-01 (34537)	—	IO 4B-A	MA-SSY48-1	—
	Standard/Standard Plus 4.8 WN (Wide Neck)	E-02 (34580)	—	IO 4B-B	MA-SSY65-1	—
	Standard Plus 3.3 NN (Narrow Neck)	G-01 (34541)	—	—***	—	—
<b>Zimmer Dental</b>	Screw-Vent / Screw-Vent friction fit 3.3, 3.7	B-01 (34527)	IO-B-01 (35212)	IO 5A-A	MA-ZSV35-1	YES
	Screw-Vent / Screw-Vent friction fit 4.5	B-02 (34528)	IO-B-02 (35213)	IO 5A-B	MA-ZSV45-1	YES
	Tapered Screw-Vent / Tapered Screw-Vent friction fit 3.5	B-01 (34527)	IO-B-01 (35212)	IO 5A-A	MA-ZSV35-1	YES
	Tapered Screw-Vent / Tapered Screw-Vent friction fit 4.5	B-02 (34528)	IO-B-02 (35213)	IO 5A-B	MA-ZSV45-1	YES
	Tapered Screw-Vent / Tapered Screw-Vent friction fit 5.7	B-03 (34529)	IO-B-03 (35214)	IO 5A-C	MA-ZSV57-1	YES
	SwissPluss 4.8	E-01 (34537)	—	—	—	—

All trademarks and company names are the property of their respective owners.

\*\*\* Scan body/ Model analog exist but is not supported for Atlantis solutions.