

Atlantis® abutments implant compatibility chart



Manufacturer	Implant	Atlantis Abutment		Atlantis Crown		Atlantis CustomBase solution - Atlantis Abutment & Atlantis Crown		Atlantis Crown Abutment		Atlantis Conus Abutment	
				Cut-back or Full-contour		Cut-back or Full-contour		Cut-back		Custom	Over-denture
		Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti*	Zirconia	Titanium Gold-shaded ti*	Angulated screw access	Titanium	Zirconia	Titanium Gold-shaded ti	Titanium
BioHorizons	INTERNAL										
	Internal 3.5, 4.5, 5.7	x	x	x	x	x	x	x	x	x	x
	Tapered Internal 3.5, 4.5, 5.7	x	x	x	x	x	x	x	x	x	x
	EXTERNAL										
BioHorizons	External implant 3.5	x								x	
	External implant 4.0, 5.0	x	x						x	x	
	INTERNAL										
	Certain 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	x	x	x	x	x		x	x	x	x
Biomet 3i	Certain Prevail 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	x	x	x	x		x	x	x	x
	Certain XP 4/5, 5/6	x	x	x	x	x		x	x	x	x
	EXTERNAL										
	Miniplant 3.25	x	x					x	x	x	x
Biomet 3i	Osseotite 3.75, 4.0, 5.0, 6.0	x	x					x	x	x	x
	Osseotite XP 4/5, 5/6	x	x					x	x	x	x
	INTERNAL										
	CONELOG 3.3	x		x		x		x		x	x
Camlog	CONELOG 3.8, 4.3	x		x		x		x		x	x
	CONELOG 5.0	x		x		x		x		x	x
	Screw-line Implant 3.3	x		x		x		x		x	x
	Screw-line / Root-line Implant 3.8, 4.3, 5.0, 6.0	x	x	x	x	x		x	x	x	x
Dentsply Sirona	INTERNAL										
	Ankylos C/X - C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)	x	x							x	
	Ankylos C/X - /X 3.5, 4.5, 5.5, 7.0 (indexed)	x	x	x	x	x	x			x	x
	OsseoSpeed EV 3.0 - Green	x		x		x	x	x			
	OsseoSpeed EV 3.6 - Purple	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed EV 4.2 - Yellow	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed EV 4.8 - Blue	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed EV 5.4 - Brown	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed Profile EV 4.2 Yellow	x	x	x		x	x	x		x	x
	OsseoSpeed Profile EV 4.8 Blue	x	x	x		x	x	x		x	x
	OsseoSpeed/OsseoSpeed TX 3.0 - Yellow	x		x		x	x	x			
	OsseoSpeed/OsseoSpeed TX 3.5, 4.0 - Aqua	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed/OsseoSpeed TX 4.5, 5.0 - Lilac	x	x	x	x	x	x	x	x	x	x
	OsseoSpeed TX Profile 4.5, 5.0	x		x		x		x		x	x
	OsseoSpeed 4.0 - Blue (Previous)	x								x	
Xive S 3.0	x		x		x	x	x				
Xive S 3.4, 3.8, 4.5, 5.5	x	x	x	x	x	x	x	x	x	x	
Hiossen	INTERNAL										
	ET Mini 3.5	x	x	x		x		x	x	x	x
Hiossen	ET Regular 4.0, 4.5, 5.0, 6.0, 7.0	x	x	x		x		x	x	x	x
	INTERNAL										
Keystone Dental	PrimaConnex SD 3.3, 3.5, Genesis 3.8	x	x	x	x	x		x	x	x	x
	PrimaConnex RD 4.0, 4.1, Genesis 4.5	x	x	x	x	x		x	x	x	x
	PrimaConnex WD 5.0, Genesis 5.5, 6.5	x	x	x	x	x		x	x	x	x

All trademarks and company names are the property of their respective owners.
 * 1 day extra turnaround time.

Please turn over →

Manufacturer	Implant	Atlantis Abutment		Atlantis Crown		Atlantis CustomBase solution - Atlantis Abutment & Atlantis Crown		Atlantis Crown Abutment		Atlantis Conus Abutment	
				Cut-back or Full-contour		Cut-back or Full-contour		Cut-back		Custom	Over-denture
		Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti*	Zirconia	Titanium Gold-shaded ti*	Angulated screw access	Titanium	Zirconia	Titanium Gold-shaded ti	Titanium
MIS	INTERNAL										
	M4 NP 3,3	x		x		x	x	x		x	x
	M4 SP 3,7, 4,2	x	x	x	x	x	x	x	x	x	x
	M4 WP 5,0, 6,0	x	x	x	x	x	x	x	x	x	x
	SEVEN NP 3,3	x		x		x	x	x		x	x
	SEVEN SP 3,75, 4,2	x	x	x	x	x	x	x	x	x	x
SEVEN WP 5,0, 6,0	x	x	x	x	x	x	x	x	x	x	
Nobel Biocare	INTERNAL										
	NobelActive 3.0	x		x		x	x	x			
	NobelActive/NobelReplace Conical Connection NP 3.5	x	x	x	x	x	x	x	x	x	x
	NobelActive/NobelReplace Conical Connection RP 4.3, 5.0	x	x	x	x	x	x	x	x	x	x
	NobelActive WP 5,5	x	x	x	x	x	x	x	x	x	x
	NobelReplace NP 3.5	x	x	x	x	x		x	x	x	x
	NobelReplace RP 4.3	x	x	x	x	x		x	x	x	x
	NobelReplace WP 5.0	x	x	x	x	x		x	x	x	x
	NobelReplace 6.0	x	x	x	x	x		x	x	x	x
	EXTERNAL										
	Brånemark System NP 3,3	x		x		x		x		x	x
	Brånemark System RP 3,75, 4,0	x	x	x		x		x	x	x	x
Brånemark System WP 5,0, 5,5	x	x	x		x		x	x	x		
Straumann	INTERNAL										
	Bone Level 3.3 NC (Narrow CrossFit)	x	x	x	x	x	x	x	x	x	x
	Bone Level 4.1, 4.8 RC (Regular CrossFit)	x	x	x	x	x	x	x	x	x	x
	Standard/Standard Plus 4.8 RN (Regular Neck)	x		x		x		x		x	x
	Standard/Standard Plus 4.8 WN (Wide Neck)	x		x		x		x		x	x
EXTERNAL											
Standard Plus 3.3 NN (Narrow Neck)	x								x		
Zimmer Dental	INTERNAL										
	Tapered Screw-Vent 3.5, 4.5, 5.7	x	x			x	x	x	x	x	x
	Tapered Screw-Vent 3.5, 4.5, 5.7 friction fit	x				x		x		x	x
	Screw-Vent 3.3, 3.7, 4.5	x	x			x	x	x	x	x	x
	Screw-Vent 3.3, 3.7, 4.5 friction fit	x				x		x		x	x
SwissPlus 4.8	x				x		x		x	x	

* 1 day extra turnaround time.

Turnaround times

Turnaround times (in-house business days) for all case sizes		
Titanium	Gold-shaded titanium	Zirconia
2	3	4

* Zirconia crowns on an Atlantis abutment in titanium or gold-shaded titanium require an extra day turnaround time.

Note: Turnaround time listed begins the first business day after receipt of all case information and material required and does not include shipping, which is determined by your choice of delivery method. Turnaround for lab-based scanned cases is one day less than the times shown.

Turnaround times are based on the following:

- Design approvals to process case must be granted by 10 am your local time.
- Lab-scanned cases, including rescans, must be uploaded to Atlantis WebOrder by 3 pm (EST) in order to save a day of processing compared with the times listed above.
- Material, on-hold and case information issues must be resolved by 8 pm (EST).

Other services

Service	Additional turnaround days
Laboratory services	add 1 day
Provisional restoration	add 5 days
Removable diagnostic preparation	add 2 days

All orders for Atlantis abutments are shipped with the corresponding abutment screw needed at no extra charge.

Note: The implant manufacturer's recommended torque should always be used when placing the abutment screw.