

DIGITAL IMPRESSION COMPARISON CHART

3M Mobile True Definition Scanner

3Shape TRIOS Scanner Wireless

CS 3600

3M ORAL CARE

3SHAPE

CARESTREAM DENTAL



SCANNER FILE OUTPUT	STL	STL, Proprietary	Proprietary, STL
SCANNER TYPE	3D Video	LED Pattern Projection	LED Pattern Projection
WAND WEIGHT	8.9 ounces	13.16 ounces	11.5 ounces
SCANNER ACCURACY (in Microns)	<5	7	11
INFECTION CONTROL	Rinse, then soak in disinfectant	Autoclavable/wipe down	Autoclavable/wipe down
TOOTH SURFACE SPRAY?	Yes	No	No
PER SCAN CLICK FEE?	No	No	No
FULLARCH SCAN TIME	1 minute	45 seconds to 2 minutes	± 1 minute
SUBSCRIPTION	\$329 monthly or \$3000/yr pay in full option	\$1800/yr starting year 2	No subscription
SCAN IN COLOR OR BW	BW	Color	Color or BW
LAPSED TIME <i>until digital model file is available</i>	10 minutes	5 minutes	5 minutes
CLOUD STORAGE	Yes	No	No
TOUCH SCREEN	Yes	Yes	No
WARRANTY?	1 year, extended warranty available	1 year, extended warranty available	3 years
IN-HOUSE TRAINING	1-2 days	1 day	1 day
SUGGESTED RETAIL PRICE	\$16,995	\$35,000 - \$70,000	\$20,000 - \$30,000
ORTHO INDICATIONS?	Yes	Yes	Yes
PARTIAL/DENTURE INDICATIONS?	No	Yes	Yes
OTHER UNIQUE FEATURES	It is the first mobile intraoral scanner that operates solely on a tablet. It is a cost effective entry into the digital dentistry.	Realistic colors and shade measurement, HD photos, documented accuracy and Patient Specific Motion	Features hybrid scan workflow, faster scanning speeds, semi-automated bite adjustment, preparation check, multiple bite capture and choice of sound options
COMPATIBLE WITH SCAN BODIES?	Yes	Yes	Yes

Prices and product descriptions in price comparisons are as published by the named suppliers on their website on the specified date. Price comparisons assume all products offered are of the same specification on the applicable date. Prices are subject to change and may vary at any time without notice. The prices shown exclude taxes, delivery charges, special offers, promotions, etc. which may alter the price payable. Price comparison information is for personal use and guidance only and does not constitute any contractual representation, warranty or obligation on our part. Liability for errors, omissions or consequential loss is expressly disclaimed all trademarks and tradenames used in price comparisons are acknowledged to be the copyright of their respective owners.

DIGITAL IMPRESSION COMPARISON CHART

CEREC AC Connect with Omnicam

iTero Element 2 Scanner

Planmeca Emerald

DENTSPLY SIRONA

ALIGN TECHNOLOGY

PLANMECA USA



	DENTSPLY SIRONA	ALIGN TECHNOLOGY	PLANMECA USA
SCANNER FILE OUTPUT	STL	STL	STL
SCANNER TYPE	Streaming	Confocal Imaging	Laser
WAND WEIGHT	11 ounces	17.6 ounces	9.1 ounces
SCANNER ACCURACY (in Microns)	25	10-25	<5
INFECTION CONTROL	Metal-Dry Heat	Disposable Plastic Sleeve - \$3	Autoclavable
TOOTH SURFACE SPRAY?	No	No	No
PER SCAN CLICK FEE?	No	Depends on Subscription	No
FULL-ARCH SCAN TIME	6-10 minute	1 minutes	45 seconds
SUBSCRIPTION	\$199 monthly, includes updates & warranty	\$360 monthly includes scans & warranty \$180 monthly includes scans Pay Per Scan Plan \$20-\$35	PlanCAD. \$4421 no annual fees
SCAN IN COLOR OR BW	Color	Color	Color
LAPSED TIME <i>until file is available</i>	5 minutes	5 minutes	0-5 minutes
CLOUD STORAGE	Yes	Yes	No
TOUCH SCREEN	No	Yes	No
WARRANTY?	1 year, extended warranty available	Depends on Subscription	3-5 years
IN-HOUSE TRAINING	1 day, additional available	Yes	Yes, depends on Distributor
SUGGESTED RETAIL PRICE	\$29,995	Varies by Region	\$30,000-\$60,000
ORTHO INDICATIONS?	Yes	Yes	Yes
PARTIAL/DENTURE INDICATIONS?	No	No	Yes
OTHER UNIQUE FEATURES	Ability to do a single-tooth restorations to full-smile rehabilitation in a single visit with continuous video option, 3D imaging and shade detection.	Uninterrupted Scanning, improved wand touchpad and HD visuals with widescreen, very popular and proven system, complete laboratory connectivity	Ergonomic design and two-button operation provide superior control and comfort, plug-and-play solution, impeccable infection control, anti-fog feature, USB-3 connection
COMPATIBLE WITH SCAN BODIES?	Yes	Yes	Yes

Prices and product descriptions in price comparisons are as published by the named suppliers on their website on the specified date. Price comparisons assume all products offered are of the same specification on the applicable date. Prices are subject to change and may vary at any time without notice. The prices shown exclude taxes, delivery charges, special offers, promotions, etc. which may alter the price payable. Price comparison information is for personal use and guidance only and does not constitute any contractual representation, warranty or obligation on our part. Liability for errors, omissions or consequential loss is expressly disclaimed all trademarks and tradenames used in price comparisons are acknowledged to be the copyright of their respective owners.