

Inspired by Nature.
Perfect by Design.

Classic★Craft
DENTAL LABORATORY

877.552.7238
classiccraftdental.com

Crowns that drive referrals

Faster Seat Times

Proven Success

Natural Shades

Easy Placement

PFZ & PFM

TWO GREAT CHOICES

Consistent

Traditional Preps

Durable

No Grey Halo at the Margin

Custom, not complicated

Perfect anatomy



Inspired by Nature.
Perfect by Design.

ZirLifeUSA™
DENTAL CROWNS GUARANTEED FOR LIFE



IPS e.max®

PORCELAIN FUSED TO ZIRCONIA

Porcelain fused to zirconia crowns are a good choice for strength, durability, and aesthetics. We take e.max powder and combine it with zirconium oxide. The crowns are layered with the e.max until the restoration is complete to your patient's case specifications. PFZ crowns are designed using our CAD/CAM software for precision and are preferred for anterior teeth. They allow more light to pass through with medium to high translucency, allowing them to be more natural compared to metal-based crowns. Our specialized intaglio groove allows for better retention after seating.

PFZ

- Up to 4 unit bridges
- Blocks out endo-treated teeth
- Medium-high translucency
- No grey halo at margin
- Intaglio groove for better retention
- Vita Classic and 3-D Master shade matching
- Conventional cement or bonding



IVOCLAR CLASSIC PORCELAIN TO ARGEN METALS

PFM

- Single unit to full round house bridges
- Porcelain is engineered to expand with the metal
- Medium translucency
- Porcelain butt margin available
- Dies are "dipped" for better retention
- Vita Classic and 3-D Master shade matching
- Conventional cement



Classic★Craft
DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007
2807-B Roger Lacey Drive, Savannah, GA 31404
classiccrafterdental.com