

Inspired by Nature.
Perfect by Design.

Classic★Craft
DENTAL LABORATORY

877.552.7238
classiccraftdental.com

Go from a slow seat to a fast seat

Put “big” back into smiles

Guaranteed for life

Drive patient referrals

Make peace
not remakes

Build a legacy practice



WITH

**Perfect and
Predictable**

ZirLifeUSA
FULL CONTOUR ZIRCONIA CROWN

YOU CAN

**Ditch that late
delivery date**

Stop the grind

**Go from frustrated
to finished early**

Don't Worry Be Happy



Stop hiding behind
“the lab usually does a better job”

Forget about failed crowns

Reclaim the control you deserve

Custom, not complicated



Inspired by Nature.
Perfect by Design.

ZirLifeUSA™
DENTAL CROWNS GUARANTEED FOR LIFE

- Indicated for posterior crowns & bridges
- >1200 mpa strength
- Medium translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 4 units
- Custom dip stained for true shades
- Polished after milling to remove micro ridges
- Lifetime warranty

ZirLifeUSA is the best solution for:

- Single crowns
- Full-arch bridges
- Screw-retained implant crowns




ZirLife™ AESTHETIC
TRANSLUCENT ZIRCONIA CROWN

- Indicated for anterior crowns & bridges
- 803 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 2 units
- Pre-shaded multi-layered zirconia for gradient of shade
- Polished after milling to remove micro ridges
- Lifetime warranty

ZirLifeAESTHETIC
is the best solution for:

- Single crowns
- Bridges with a maximum of two pontics
- In-lays
- On-lays



Classic★Craft
DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007
2807-B Roger Lacey Drive, Savannah, GA 31404
classiccraftdental.com