

Implant Connection Guide

Implant System	Generation 1	Generation 2	Available Sizes				
ASTRA TECH™ EV*	SEV Series	-	3.0	3.6	4.2	4.8	5.4
ASTRA TECH™ OsseoSpeed TX*	S Series	AO Series	3.5/4.0		4.5/5.0		
BioHorizons®*	BIO Series	-	3.0	3.5	4.5	5.7	
Biomet 3i™ Certain*	H Series	BC Series	3.4 NP	4.1 RP	5.0 WP		
Biomet 3i™ OSSEOTITE*	I Series	BO Series	3.5	4.1	5.0		
CAMLOG®*	CAM Series	-	3.3	3.8	4.3	5.0	
Dentsply Friadent XiVE®*	T Series	DX Series	3.4	3.8	4.5	5.5	
Hiossen ET / HG®*	HIO Series	-	3.5	4.0/4.5/5.0/6.0			
Neoss®*	NEO Series	-	3.5/4.0/4.5/5.0/5.5				
NobelActive™*	F Series	NC Series	3.5 NP	4.3/5.0 RP			
NobelReplace Select®*	E Series	NR Series	3.5 NP	4.3 RP	5.0 WP	6.0	
Nobel Brånemark®*	K Series	NB Series	3.5	4.1	5.1		
Straumann Bone Level®*	L Series	SB Series	3.3 NC	4.1/4.8 RC			
Straumann SynOcta®*	N Series	SS Series	3.5 NN	4.8 RN	6.5 WN		
Zimmer Tapered Screw-Vent®*	R Series	ZS Series	3.5	4.5	5.7		

* Not trademarks of the Argen Corporation.

