

What does this mean for our customers?

### More options, including more screw-retained capabilities


- Titanium bases can be customized 4-10mm tall so we can maximize the crown to Ti-base ratio
- Engaging and non-engaging options for single units and bridges
- Choose screw-retained or custom hybrid zirconia restorations
- No intraoral cementing

### Save time and money

- Faster turnaround time versus custom titanium abutments
- We do these restorations in-house
- More economical

### Usability & accurate scanning

- Screw head is compatible with the implant manufacturer driver
- Scan bodies contain a screw that will not fall out of the scan body in a patient's mouth
- The scan body height is comfortable for the patient
- The interface is metal so that you can confirm that it's properly seated with a radiograph
- The coronal portion is PEEK, which is excellent for scanning accuracy



### Why does Classic Craft Dental Laboratory choose Open Implants?

Open Implants is a **US manufacturer** of prosthetic components who has designed a system for the fabrication of esthetic restorations in today's digital world. OI has **FDA 510K clearance** on their components, which are **compatible with the top implant systems**. They have **high quality standards** and they stand by their quality with a **warranty on their products as well as the implant**.

### Compatible Brands:

AstraTech® Osseospeed  
Biomet® 3i Certain  
Nobel Biocare® Replace Select  
Nobel Biocare® Nobel Active  
Nobel Biocare® Multi Unit  
Straumann® Bone Level  
Zimmer® Tapered Screw Vent

### For more information, contact:

Trevor Gottfredson  
CCDL Implant Specialist  
877.552.7238 Office

You can also visit  
[www.openimplants.com](http://www.openimplants.com)  
or call 1-781-587-3242

