


# IMPLANT BARS



Classic Craft Dental Lab now offers milled titanium or custom-cast bars for over dentures. We have partnered with Atlantis and Panthera to offer a unique variety of options for cases that are more complex in nature. Custom bars can be paired with a semi-fixed or fixed prosthesis as dictated by your treatment plan. Custom-cast locator and Hader bars are also available.

## CUSTOM-CAST BAR FEATURES:

- Solid titanium for increased strength
- Increased flexibility when addressing A-P spread issues
- Multiple attachment options
- Can be cast with engaging or non-engaging UCLA abutments

## CUSTOM-CAST BARS ARE THE BEST SOLUTION FOR:

- Increasing load distribution over multiple implants
- Correcting difficult implant angles
- Patient treatment plans requiring semi-fixed prosthesis



**DIAMART**  
*Domenico Cascione*  
BY  
DOMENICO CASCIONE, CDT



**Atlantis**<sup>®</sup>  
Provided by Dentsply Sirona Implants

If you have questions call Trevor Gottfredson at  
Classic Craft Dental Laboratory at 877.552.7238 for assistance.