

Inspired by Nature.  
Perfect by Design.

**Classic★Craft**  
DENTAL LABORATORY

877.552.7238  
classiccraftdental.com

**Go from a slow seat to a fast seat**

**Put "big" back into smiles**

Guaranteed for life

**Drive patient referrals**

Make peace  
not remakes



Build a legacy practice

**WITH**

**Perfect and  
Predictable**

**ZirLifeUSA**  
FULL CONTOUR ZIRCONIA CROWN

**YOU CAN**

**Ditch that late  
delivery date**

**Stop the grind**

**Go from frustrated  
to finished early**

**Don't Worry Be Happy**



Stop hiding behind  
"the lab usually does a better job"

Forget about failed crowns

Reclaim the control you deserve

**Custom, not complicated**




Inspired by Nature.  
Perfect by Design.

**ZirLifeUSA™**  
DENTAL CROWNS GUARANTEED FOR LIFE

- Indicated for posterior crowns & bridges
- >1200 mpa strength
- Medium translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 4 units
- Custom dip stained for true shades
- Polished after milling to remove micro ridges
- Lifetime warranty

**ZIRLIFE USA** is the best solution for:

- Single crowns
- Full-arch bridges
- Screw-retained implant crowns



**Classic★Craft**  
DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007

2807-B Roger Lacey Drive, Savannah, GA 31404

classiccrafterdental.com

©2020 CLASSIC CRAFT DENTAL LABORATORY

**ZirLife™ AESTHETIC**  
TRANSLUCENT ZIRCONIA CROWN

- Indicated for anterior crowns & bridges
- 803 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 2 units
- Pre-shaded multi-layered zirconia for gradient of shade
- Polished after milling to remove micro ridges
- Lifetime warranty

**ZIRLIFE AESTHETIC**

is the best solution for:

- Single crowns
- Bridges with a maximum of two pontics
- In-lays
- On-lays

