

Inspired by Nature.
Perfect by Design.

ZirLifeUSA™

FULL CONTOUR ZIRCONIA CROWN

CHAIRSIDE TABLE GUIDE

- Indicated for posterior crowns & bridges
- >1100 mpa strength
- Medium Translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Dipped stained for true shades
- Polished after milling to remove micro ridges




©2020 CLASSIC CRAFT DENTAL LABORATORY

Classic★Craft

DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007
2807-B Roger Lacey Drive, Savannah, GA 31404
classiccraftdental.com



POSTERIOR CROWNS

SEE ANTERIOR PREP ON OTHER SIDE




Inspired by Nature.
Perfect by Design.

ZirLife™ AESTHETIC


TRANSLUCENT ZIRCONIA CROWN

CHAIRSIDE TABLE GUIDE

- Indicated for anterior crowns & bridges
- 803 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Dipped stained for true shades
- Polished after milling to remove micro ridges



ANTERIOR CROWNS



ANTERIOR BRIDGES

SEE POSTERIOR PREP ON OTHER SIDE