

Inspired by Nature.  
Perfect by Design.

# IMPLANT BARS

Classic Craft Dental Lab now offers milled titanium or custom-cast bars for over dentures. We have partnered with Atlantis and Panthera to offer a unique variety of options for cases that are more complex in nature. Custom bars can be paired with a semi-fixed or fixed prosthesis as dictated by your treatment plan. Custom-cast locator and Hader bars are also available.



Atlantis®

Provided by Dentsply Sirona Implants



## CUSTOM-CAST BAR FEATURES:

- Solid titanium for increased strength
- Increased flexibility when addressing A-P spread issues
- Multiple attachment options
- Can be cast with engaging or non-engaging UCLA abutments

### CUSTOM-CAST BARS

#### ARE THE BEST SOLUTION FOR:

- Increasing load distribution over multiple implants
- Correcting difficult implant angles
- Patient treatment plans requiring a semi-fixed prosthesis

If you have questions call Trevor Gottfredson  
at Classic Craft Dental Laboratory  
at 877.552.7238 for assistance.

**Classic★Craft**  
DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007

2807 Roger Lacey Drive Ste. B, Savannah, GA 31404

classiccraftdental.com

©2021 CLASSIC CRAFT DENTAL LABORATORY



DIAMART™  
Domenico Cascione  
BY  
DOMENICO CASCIONE, CDT

Patent pending

