

Inspired by Nature.
Perfect by Design.

ZirLifeUSA™

FULL CONTOUR ZIRCONIA CROWN

CHAIRSIDE TABLE GUIDE

- Indicated for posterior crowns & bridges
- >1200 mpa strength
- Medium Translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 4 units
- Pre-shaded for predictable, consistent shades
- Polished after milling to remove micro ridges
- Lifetime warranty




©2021 CLASSIC CRAFT DENTAL LABORATORY

Classic★Craft

DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007
2807 Roger Lacey Drive Ste. B, Savannah, GA 31404
classiccraftdental.com



POSTERIOR CROWNS

SEE ANTERIOR PREP ON OTHER SIDE



Inspired by Nature.
Perfect by Design.

ZirLife™ AESTHETIC

TRANSLUCENT ZIRCONIA CROWN

CHAIRSIDE TABLE GUIDE


- Indicated for anterior crowns & bridges
- 900 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 2 units
- Pre-shaded multi-layered zirconia for gradient of shade
- Polished after milling to remove micro ridges
- Lifetime warranty




©2021 CLASSIC CRAFT DENTAL LABORATORY

Classic★Craft
DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007
2807 Roger Lacey Drive Ste. B, Savannah, GA 31404
classiccraftdental.com



ANTERIOR CROWNS



ANTERIOR BRIDGES

SEE POSTERIOR PREP ON OTHER SIDE

