




Image courtesy of Dentsply Sirona, Inc.

FREE SCAN BODIES

Classic Craft Dental Lab has partnered with Atlantis to provide free scan bodies for a wide selection of implant manufacturers. Scan bodies can be used on all digital platforms. Atlantis scan bodies create a seamless transfer of unique implant placement information to us for an unmatched custom abutment and aesthetic crown result.



Atlantis®
Provided by Dentsply Sirona Implants

ABUTMENTS



Image courtesy of Dentsply Sirona

CUSTOM ABUTMENTS

Classic Craft Dental Lab works with all custom abutment manufacturers including, but not limited to: Atlantis, Argen, Zimmer Biomet, Straumann, Nobel Biocare, and Vulcan. Together we produce abutments and implant crowns that meet or exceed the specifications of your treatment plan. We are also an Encode Empowered lab, which allows us to unlock Encode healing caps for abutment milling without using a third party.

CUSTOM ABUTMENT FEATURES:

- Custom abutment partners assume all warranties for OEM implants
- Options include: titanium, gold hue titanium, zirconia, hybrid zirconia, and one piece, screw-retained on a ti-base

CUSTOM ABUTMENTS ARE THE BEST SOLUTION FOR:

- Difficult implant placement
- Enhanced tissue management
- Natural contours
- Screw-retained crowns
- Angle correction

OEM/STOCK ABUTMENTS ARE THE BEST SOLUTION FOR:

- Implants that are packaged with an abutment for cost savings
- Lesser known implant companies
- Restorations that don't require special needs for tooth or tissue management

OEM - STOCK

(ORIGINAL EQUIPMENT MANUFACTURERS)

Classic Craft Dental Lab works with all FDA-approved implant manufacturers and each case we produce with an OEM or stock abutment comes with a soft tissue model and placement jig. If minor modifications are needed for better tooth and tissue contours, they are included in the implant fee. We will never substitute an OEM abutment without your approval and will always work with you to determine the best option for restoring the implant.



877.552.7238 Fax: 912.355.1007
classiccraftdental.com

©2022 CLASSIC CRAFT DENTAL LABORATORY



IMPLANT BARS

Classic Craft Dental Lab now offers milled titanium or custom-cast bars for over dentures. We have partnered with Atlantis and Panthera to offer a unique variety of options for cases that are more complex in nature. Custom bars can be paired with a semi-fixed or fixed prosthesis as dictated by your treatment plan. Custom-cast locator and Hader bars are also available.



Atlantis®

Provided by Dentsply Sirona Implants



CUSTOM-CAST BAR FEATURES:

- Solid titanium for increased strength
- Increased flexibility when addressing A-P spread issues
- Multiple attachment options
- Can be cast with engaging or non-engaging UCLA abutments

CUSTOM-CAST BARS ARE THE BEST SOLUTION FOR:

- Increasing load distribution over multiple implants
- Correcting difficult implant angles
- Patient treatment plans requiring a semi-fixed prosthesis



If you have questions call our Implant Specialist at Classic Craft Dental Laboratory at 877.552.7238 for assistance



877.552.7238 Fax: 912.355.1007
classiccraftdental.com

©2022 CLASSIC CRAFT DENTAL LABORATORY

