



# REMOVABLES

As the demand for expert removable care continues to increase, so does the importance of partnering with a quality lab. Getting a denture can restore a smile and change a patient's way of life. We offer standard and premium full dentures, frameworks, over-denture milled bars, All on 4 restorations. We also offer a full line of removable services, such as denture repairs, adding teeth, and hard relines.

- Premium full dentures
- Standard full dentures
- Frameworks
- Over-denture milled bars
- All on 4 restorations





# DURACETAL™

COSMETIC APPLIANCES

DuraCetal® is a pure, highly crystalline acetal copolymer resin possessing high tensile and flexural strength, fatigue resistance, and hardness. These properties, when combined with its very low moisture absorption, make Myerson DuraCetal® an ideal, light weight alternative to cast metal partial denture frameworks. It is also indicated for clasps and unilateral partials. For patients that do not want a metal framework, DuraCetal® clasps and frameworks are the solution.

## BENEFITS OF DURACETAL®

- More rigid, flexible & comfortable than traditional metal frameworks
- Bacteria and stain resistant with no odor
- Better fit due to flex from "snap fit" feature
- Natural colors: Opaque Tooth Colored or Clear
- Snap-fit with excellent retention.
- No adjusting as clasps flex & return to original shape.
- Biocompatible, monomer free, and FDA approved
- Strong and durable
- Can be "added to" or repaired
- Can be adjusted with a carbide bur
- Excellent esthetics – looks natural in the mouth



## ~ HYBRID PARTIALS

- Aesthetic clasp
- Metal framework on the lingual for stability
- Teeth can be added
- Available with acrylic saddles

## VALPLAST PARTIALS ›

- Comfortable and light
- Aesthetic - no metal show through
- Longevity - very durable



## ‹ CAST PARTIALS

- All frameworks are Ni & Be free
- Metal framework on the lingual for stability
- Teeth can be added
- Acrylic saddle(s)
- Very durable

