

Go from a slow seat to a fast seat



Guaranteed for life

Put "big" back into smiles

Drive patient referrals

Make peace
not remakes

Build a legacy practice



Stop the grind

WITH

Perfect and
Predictable

ZirLifeUSA
FULL CONTOUR ZIRCONIA CROWN

YOU CAN

Ditch that late
delivery date

Don't Worry Be Happy


Go from frustrated
to finished early

Forget about failed
crowns

Stop hiding behind
"the lab usually does a better job"

Reclaim the control you deserve

Custom, not complicated



- Indicated for posterior crowns & bridges
- >1200 mpa strength
- Medium translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 4 units
- Pre-shaded for predictable, consistent shades
- Polished after milling to remove micro ridges
- Lifetime warranty

ZIRLIFE USA is the best solution for:

- Single crowns
- Full-arch bridges
- Screw-retained implant crowns



*Inspired by Nature.
Perfect by Design.*



ZirLife™ AESTHETIC
TRANSLUCENT ZIRCONIA CROWN

- Indicated for anterior crowns & bridges
- 900 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 2 units
- Pre-shaded multi-layered zirconia for gradient of shade
- Polished after milling to remove micro ridges
- Lifetime warranty

ZIRLIFE AESTHETIC

is the best solution for:

- Single crowns
- Bridges with a maximum of two pontics
- In-lays
- On-lays

