


## CHAIRSIDE TABLE GUIDE



POSTERIOR CROWNS



- Indicated for posterior crowns & bridges
- >1200 mpa strength
- Medium Translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 4 units
- Pre-shaded for predictable, consistent shades
- Polished after milling to remove micro ridges
- Lifetime warranty



**Classic★Craft**  
DENTAL LABORATORY


877.552.7238 Fax: 912.355.1007  
classiccraftdental.com

©2022 CLASSIC CRAFT DENTAL LABORATORY




## CHAIRSIDE TABLE GUIDE

- Indicated for anterior crowns & bridges
- 900 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 2 units
- Pre-shaded multi-layered zirconia for gradient of shade
- Polished after milling to remove micro ridges
- Lifetime warranty



ANTERIOR CROWNS



ANTERIOR BRIDGES



**Classic★Craft**  
DENTAL LABORATORY

877.552.7238 Fax: 912.355.1007

[classiccrafterdental.com](http://classiccrafterdental.com)

©2022 CLASSIC CRAFT DENTAL LABORATORY

