

## ZirLife™

## ZirLifeUSA™

DENTAL CROWNS GUARANTEED FOR LIFE



- >1200 mpa strength
- Medium translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 4 units
- Polished after milling to remove micro ridges
- Lifetime warranty

## ZirLife™ AESTHETIC

TRANSLUCENT ZIRCONIA CROWN

- 900 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 2 units
- Polished after milling to remove micro ridges
- Lifetime warranty



## IPS e.max®

### PRESSED ANTERIOR

- 450 mpa strength
- 97% clinical rating
- 96.8% survival rate

### MILLED POSTERIOR

- 500 mpa strength
- 97% clinical rating
- 96.8% survival rate

IVOCLAR  
CLASSIC PORCELAIN  
TO ARGEN  
METALS (PFM)

ZirLifeUSA™  
+ IPS e.max®

PORCELAIN FUSED  
TO ZIRCONIA (PFZ)

FULL - CAST  
RESTORATIONS



## Classic★Craft

DENTAL LABORATORY