

CLEAR COMFORT BITE SPLINT ›

- Hard night guard fabricated from modified composite material.
 - Placed in warm water before inserting in the mouth.
 - Does not contain methyl methacrylate.
- Chairside repairs can easily be made with cold-cure acrylics.



CLASSIC COMFORT H/S™ › BITE SPLINT

- Made of two layers: a 1mm soft polyurethane inner layer and a 3mm hard copolyester outer layer.
- Standard feature with a flat occlusal plane and slight opposing cusp indentation or choose your specific bite plan needs.
 - Fabricated from ISO-certified ERKODENT® copolyester and polyurethane discs.
 - BPA-free and non-allergenic.




HARD ACRYLIC BITE SPLINT

- For patients who grind and may be showing early signs of wear.
 - Vacuum formed Erkoform material.
 - Lightweight and clear.

Bite Splint Care: Brush with a clean non-abrasive gel toothpaste.

Argen
Clear**Aligners**



GET STARTED

Scan to visit CCDL
Argen Clear Aligners.



Classic Craft
DENTAL LABORATORY