

What does this mean for our customers?

More options, including more screw-retained capabilities


- Titanium bases can be customized 4-10mm tall so we can maximize the crown to Ti-base ratio
- Engaging and non-engaging options for single units and bridges
- Choose screw-retained or custom hybrid zirconia restorations
- No intraoral cementing

Save time and money

- Faster turnaround time versus custom titanium abutments
- We do these restorations in-house
- More economical

Usability & accurate scanning

- Screw head is compatible with the implant manufacturer driver
- Scan bodies contain a screw that will not fall out of the scan body in a patient's mouth
- The scan body height is comfortable for the patient
- The interface is metal so that you can confirm that it's properly seated with a radiograph
- The coronal portion is PEEK, which is excellent for scanning accuracy



Why does Classic Craft Dental Laboratory choose Open Implants?

Open Implants is a **US manufacturer** of prosthetic components who has designed a system for the fabrication of esthetic restorations in today's digital world. OI has **FDA 510K clearance** on their components, which are **compatible with the top implant systems**. They have **high quality standards** and they stand by their quality with a **warranty on their products as well as the implant**.

Compatible Brands:

AstraTech® Osseospeed
Biomet® 3i Certain
Nobel Biocare® Replace Select
Nobel Biocare® Nobel Active
Nobel Biocare® Multi Unit
Straumann® Bone Level
Zimmer® Tapered Screw Vent

For more information, contact:

Tyler Arnold
CCDL Implant Specialist
tyler@classiccrafterdental.com
877.552.7238 Laboratory

You can also visit
www.openimplants.com
or call 1-781-587-3242

