

ZirLife™

CROWN & BRIDGE

ZirLifeUSA™

DENTAL CROWNS GUARANTEED FOR LIFE




- >1200 mpa strength
- Medium translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 4 units
- Polished after milling to remove micro ridges
- Lifetime warranty

ZirLife™ AESTHETIC

TRANSLUCENT ZIRCONIA CROWN

- 900 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 2 units
- Polished after milling to remove micro ridges
- Lifetime warranty



IPS e.max®

PRESSED ANTERIOR

- 450 mpa strength
- 97% clinical rating
- 96.8% survival rate

MILLED POSTERIOR

- 500 mpa strength
- 97% clinical rating
- 96.8% survival rate

IVOCLAR
CLASSIC PORCELAIN
TO ARGEN
METALS (PFM)

ZirLifeUSA™
+ IPS e.max®
PORCELAIN FUSED
TO ZIRCONIA (PFZ)

FULL - CAST
RESTORATIONS

 **Classic★Craft**
DENTAL LABORATORY

877.552.7238 | pics@classiccrafterdental.com | classicrafterdental.com | 2807 Roger Lacey Drive Ste. B, Savannah, GA 31404