

~ CLEAR COMFORT BITE SPLINT

- Hard night guard fabricated from modified composite material.
- Moldable after placing in warm water.
- Does not contain methyl methacrylate.
- Chairside repairs can easily be made with cold-cure acrylics.

CLASSIC COMFORT H/S™ ~ BITE SPLINT

- Made of two layers: a 1mm soft polyurethane inner layer and a 3mm hard copolyester outer layer.
- Comes Standard with a flat occlusal plane and slight opposing cusp indentation or choose your specific bite plan needs.
 - Fabricated from ISO-certified ERKODENT® copolyester and polyurethane discs.
 - BPA-free and non-allergenic.



OTHER OFFERINGS:

Hawley Retainers
Spring Retainers
Lower Lingual Arch
Band & Loop
Essix Retainer
Sports Guard
Bleaching Tray Set

HARD ACRYLIC BITE SPLINT ~

- For patients who grind and may be showing early signs of wear.
- Vacuum formed Erkoform material.
- Lightweight and clear.

WHAT TO EXPECT

Now more than ever, patients are asking for a clear aligner solution. With these aligners you can expect:


- Reliable customer service and support for every case, from file upload to aligner delivery.
- Average refinement rate of 5% and up to 20% fewer aligner stages results in an ideal treatment plan.
- Intuitive software that is easy for experienced and new staff to navigate.
- Made in the USA with premium thermoformed materials.
- Three options: Xpress (12 Clear Aligners), Limited (24 Clear Aligners) or Unlimited (25+ Clear Aligners) with complimentary retainer included in all cases.



SMILE DESIGN STUDIO

We utilize real-time automatic adjustments and comprehensive tooth-movement algorithms to produce ultimate treatment outcomes.

The intuitive user interface keeps control in your hands and offers the flexibility to fully customize every patient's case based on your treatment preferences. Attachments, interproximal reduction (IPR), final archform shape, and tooth movement and rotation speed can all be tailored to fit your standards.



\$25
OFF

Your first
Aligner Treatment

USE CODE
"ALIGNER25"
IN NOTES SECTION

GET STARTED

Scan to visit Classic Craft
Clear Aligners and
SmileDesign Studio login



Powered by G&H Orthodontics®