




DURACETAL™

COSMETIC APPLIANCES

DuraCetal® is a pure, highly crystalline acetal copolymer resin possessing high tensile and flexural strength, fatigue resistance, and hardness. These properties, when combined with its very low moisture absorption, make Myerson DuraCetal® an ideal, light weight alternative to cast metal partial denture frameworks. It is also indicated for clasps and unilaterals partials. For patients that do not want a metal framework, DurAcetal® clasps and frameworks are the solution.

BENEFITS OF DURACETAL®:

- ❖ More rigid, flexible & comfortable than traditional metal frameworks
- ❖ Bacteria and stain resistant with no odor
- ❖ Better fit due to flex from "snap fit" feature
- ❖ Natural colors: Clear or Opaque Tooth Colored
- ❖ Snap-fit with excellent retention. No adjusting as clasps flex & return to original shape.
- ❖ Bio compatible, monomer free, and FDA approved
- ❖ Strong and durable
- ❖ Can be "added to" or repaired
- ❖ Can be adjusted with a carbide bur
- ❖ Excellent esthetics – looks natural in the mouth



Classic★Craft
DENTAL LABORATORY