


OTHER OFFERINGS: Hawley Retainer, Spring Retainer, Sports Guards, Lower Lingual Arch, Band & Loop, Essix Retainer, and Bleaching Tray Set
Bite Splint Care: Brush with a clean, non-abrasive gel toothpaste.



CLEAR COMFORT BITE SPLINT ~

- ❖ Hard night guard fabricated from modified composite material.
- ❖ Moldable after placing in warm water.
- ❖ Does not contain methyl methacrylate.
- ❖ Chairside repairs can easily be made with cold-cure acrylics.



~ CLASSIC COMFORT H/S™ BITE SPLINT

- ❖ Made of two layers: a 1mm soft polyurethane inner layer and a 3mm hard copolyester outer layer.
- ❖ Comes Standard with a flat occlusal plane and slight opposing cusp indentation or choose your specific bite plan needs.
- ❖ Fabricated from ISO-certified ERKODENT® copolyester and polyurethane discs.
- ❖ BPA-free and non-allergenic.



~ HARD ACRYLIC BITE SPLINT

- ❖ For patients who grind and may be showing early signs of wear.
- ❖ Vacuum formed Erkoform material.
- ❖ Lightweight and clear.




EMA

- ❖ Smallest of our sleep appliances.
- ❖ Opens the bite while advancing the mandible to increase the airway.
- ❖ 2mm vacuum form material joined with straps of varying lengths and flexibility.
- ❖ Free of metals and the elastic straps can be changed

DREAMTAP

- ❖ Fully adjustable on a horizontal plane.
- ❖ Separate upper and lower trays.
- ❖ Adjustable screw and hook which interfaces with a socket or bar to hold the mandible in the desired position.





CLASSIC


SMILE ALIGNERS

WHY OUR CLEAR ALIGNERS?

- ❖ Made from a premium thermoformed plastic material that is engineered for precision fit, ongoing force retention for faster treatment, exceptional durability, and greater comfort for your patient.
- ❖ Tear and stain resistant.
- ❖ Manufactured in the USA using FDA 510(k) cleared, MDSAP certified processes.
- ❖ Each clear aligner is marked with a unique identifier code for the case and the stage.



STREAMLINED FOR EASY CASE MANAGEMENT AND ORDER



DOWNLOAD OUR FORMS

SCAN HERE, VISIT OUR WEBSITE, OR CALL THE LAB



CLASSIC SMILE ALIGNER PACKAGES

Xpress (12 Clear Aligners) - 1 revision, attachments & retainer included in package	\$795
Limited (24 Clear Aligners) - 1 revision, attachments & retainer included in package	\$1,095
Unlimited (25+ Clear Aligners) - 1 revision, attachments & retainer included in package	\$1,495
Clear Aligners Design Only	\$289
Per Aligner Package (Minimum of 5 Clear Aligners) - Retainer not included	\$70
Replacement Clear Aligner	\$70



Classic Craft
DENTAL LABORATORY

Powered by Argen®