

Premium full dentures

Digital dentures

Standard full dentures

Copy full dentures


Frameworks

Over-denture milled bars

All-on-4 restorations

Partials

REMOVABLES



DIGITAL DENTURES

Digital dentures are a new way to deliver predictable fit and comfort to your patients. These dentures will cut down on chair time by following a 3-step protocol. The digital design process allows for changes to be made by the click of a mouse. The material has a brilliant finish for maximum aesthetics while offering superior strength to processed dentures. After your patient receives the digital denture, a file is stored to be able to duplicate this denture in the future for peace of mind or as a back up for your patient to keep. There are multiple processes available to achieve the finished product, whether your patient needs their first denture, an exact copy of their existing denture, or an improvement on their existing denture. Start your journey into digital dentures by contacting your sales rep or our CAD/CAM department today!



CAST PARTIALS

All frameworks are Ni & Be free
Acrylic saddle(s)



VALPLAST PARTIALS

Tissue colored framework



HYBRID PARTIALS

Metal framework on the lingual for stability
Available with acrylic saddles



BENEFITS

- ❖ CAD/CAM can produce precise fit within 100 microns.
- ❖ Fewer try-in appointments. Final denture delivered in 3 appointments or less.
- ❖ Fewer post-insertion adjustment visits.
- ❖ Digital record makes it easy to replace lost or damaged dentures with a copy denture.
- ❖ Two times stronger than traditional acrylic.
- ❖ Can be repaired with acrylic.
- ❖ Made with premium Lucitone Digital Print Denture System by Dentsply Sirona.

DURACETAL™
COSMETIC APPLIANCES



DURACETAL BENEFITS

- ❖ Opaque Tooth Colored or Clear
- ❖ Can be "added to" or repaired
- ❖ Snap-fit with excellent retention
- ❖ Bacteria and stain resistant
- ❖ Strong and durable
- ❖ Biocompatible, monomer free, and FDA approved



Classic★Craft
DENTAL LABORATORY



877.552.7238



pics@classiccraftdental.com



classiccraftdental.com



2807 Roger Lacey Drive Ste. B, Savannah, GA 31404